**Document downloaded from:** <https://www.mainstaymedical.com>

**Document Information:**

Device Name: Mainstay Tunneler: Model TUN1

Document/Label Type: Instructions For Use

CE/ UK/ Australia


Language: Multi Language

Revision: C

Filename: Mainstay Tunneler—IFU—CE/UK/ Australia—Multi Language—Rev C.pdf

## **Mainstay Tunneler: Model TUN1** Instructions For Use

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## 1. DESCRIPTION

The Mainstay Tunneler Model TUN1 is designed for use during implantation of Mainstay Medical leads. It is recommended for subcutaneous tunneling between incisions in the back.

The Tunneler consists of five components: a stainless steel shaft with a formed handle, two Pebax (7233 Clear) sleeves and two stainless steel tips.

## 2. SYMBOLS AND DEFINITIONS

Do not re-sterilize

Do not re-use

Do not use if package is damaged

Consult instructions for use

Date of manufacture

Caution

Sterilized using ethylene oxide

Use-by date

Medical device

Single sterile barrier system with protective packaging inside

Unique device identifier

Contains Hazardous Substances

Catalogue number

Outside diameter

Length

Lot number

Quantity

Keep dry

Temperature limit

CE mark

UKCA mark

### 3. INTENDED PURPOSE / INDICATIONS FOR USE

The Tunneler is intended for use to route lead(s) between incisions in the back.

Tunneler is intended to be used by an implanting clinician during ReActiv8 IPG/ Leads implant surgery to make a pathway through a patient's tissue to allow the ReActiv8 percutaneous stimulation leads to be connected to the ReActiv8 IPG.

### 4. INTENDED USERS / TARGET POPULATION

The Mainstay Tunneler is intended to be used by an implanting clinician during ReActiv8 IPG/Leads implant surgery to make a pathway through a patient's tissue to allow the ReActiv8 Percutaneous Stimulation Leads to be connected to the ReActiv8 IPG.

### 5. CONTRAINDICATIONS, WARNINGS

There are no contraindications or warnings related to this device.

### 6. PRECAUTIONS

The Tunneler is a single-use-only device. Never resterilize or reuse it.

Do not resterilize the Tunneler because functionality and sterility cannot be ensured, and infections may occur. Return any open packages that have not been used to the manufacturer.

Extra Tunnelers should be available in case sterility is compromised or if the Tunneler is damaged or dropped during surgery.

Care should be taken not to injure any nerves, arteries, or veins during the tunneling procedure.

## 7. PRODUCT MATERIALS

The Tunneler is manufactured with materials that come in contact with tissue. The following materials are intended to come in contact with tissue.

- 304L Stainless Steel
- PEBAX 7233 SA 01 MED Natural


This device contains the following substance(s) defined as CMR 1B in a concentration above 0.1% weight by weight: Cobalt; CAS No. 7440-48-4; EC No. 231-158-0.

## 8. STORAGE AND HANDLING

Do not store the Tunneler where it is exposed to water or other liquids, because moisture can damage the seal integrity of the package materials.

Prior to opening the sterile package, examine it carefully for damage or compromised sterility. If the outer package has been opened or damaged, sterility of the Tunneler cannot be guaranteed and should not be used. If the package has evidence of exposure to extreme temperatures or moisture, or if there is any indication of external damage, the package should be left unopened and returned to the manufacturer.

Store the Tunneler between -29°C and 60°C.

To open the Tunneler package, do the following:

1. Carefully break tamper evident sticker on sales package and remove the sterile packaged Tunneler. The sterile package is configured as a tray in pouch.
2. Open the outer pouch.
3. Using sterile technique, remove the sterile inner tray.
4. Grasp the inner tray's tab and carefully peel off the cover to expose the contents without dropping them.

Make sure to remove all five pieces from the tray (1 shaft, 2 sleeves, and 2 tips).

## 9. DIRECTIONS FOR USE

The following directions cover the use of the Tunneler.

Do not use the Tunneler if the following occurs:

- The outer or inner package has been pierced or altered, because this can render it nonsterile.
- The expiration (use-by) data has passed, because this can adversely

- affect the component's sterility.
- Do not resterilize the Tunneler because functionality and sterility cannot be ensured, and infections may occur. Return any open packages that have not been used to the manufacturer.

### Assemble the Tunneler

The Tunneler must be assembled prior to being used. To do so, follow the steps below:

1. Using sterile technique, remove the five Tunneler components from the inner tray and place them into the sterile field.
2. Slide one of the sleeves over the shaft.
3. Carefully screw the appropriate tip onto the shaft.

Do not over tighten the tip. Doing so could make it difficult to remove.

### 9.2 Using the Tunneler

After the incisions have been made and the lead(s) have been implanted, the Tunneler can be inserted and passed from one incision to the other. If necessary, the Tunneler can be manually shaped to help direct it through the back.

Do not shape the Tunneler more than 25 degrees, because it may cause the sleeve to kink or become difficult to remove.

To pass the Tunneler, do the following:

1. Place the tip end of the Tunneler through an incision exerting force on the handle and directing the Tunneler as necessary.

Take care if tunneling towards an incision near lead; impact between the Tunneler and the lead may cause damage to the lead.

2. After the tip has passed from an incision to the other, unscrew the tip and withdraw the shaft, leaving the sleeve extending through both incisions.
3. Carefully insert the lead(s) into the sleeve until visible at the opposite incision.
4. Carefully remove the sleeve from the lead(s) at the exiting incision.
5. When the Tunneler is no longer needed, discard the Tunneler.

Used Tunneler components may be contaminated with potentially infectious substances of human origin and are considered a sharp. Utilize proper biohazard and sharps disposal.

## 10. INFORMATION AND SUPPORT

If there are questions regarding the use of the Tunneler, please contact:

A copy of the Summary of Safety and Clinical Performance (SSCP) can be viewed by searching the device brand name on the Eudamed website at <https://ec.europa.eu/tools/eudamed> when it is available. Electronic Instructions for Use (eIFU) including a copy of the SSCP is also available at [www.mainstaymedical.com/resources](http://www.mainstaymedical.com/resources).

Please report any serious incident associated with the Mainstay product to your national competent authority and to Mainstay at the following:

Via postal mail:

Mainstay Medical Limited  
Clonmel House  
Forster Way  
Swords, K67F2K3 Co. Dublin  
Ireland

Internet: [www.mainstaymedical.com](http://www.mainstaymedical.com)

Via email: [contact@mainstaymedical.com](mailto:contact@mainstaymedical.com)

Via phone: +353 766 801428

UK Responsible Person:

ICON (LR) Limited  
500 South Oak Way  
Green Park, Reading RG2 6AD  
United Kingdom

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Mainstay Medical (Australia) Pty Ltd