Document downloaded from: <https://www.mainstaymedical.com>

Document Information

Document Number: 990082-001

Device Name: ReActiv8 Release Tool

Document/Label Type: Instructions for Use

Country: Global

Language: Multilanguage

Revision: C

Filename: 990082-001_ReActiv8 Release Tool IFU (Global) - Multilanguage_Rev C

(English)

ReActiv8 Release Tool: Model EXT1

Intended Purpose/Indications for Use

The ReActiv8 Release Tool (REF EXT1) is designed for use during the removal of the ReActiv8 leads.

Clinical Benefits

Patients may benefit by way of reduced risk of unsuccessful lead removal and abandonment of the indwelling lead or lead components.

Warnings & Precautions

- Do not use excessive force as this can cause the lead to fracture, potentially resulting in lead fragmentation.
- Do not cut the lead too short as there will be insufficient length to facilitate use of the tool.
- When disconnecting the lead from the IPG, take care to leave sufficient length of the lead so that the proximal end can still be grasped when the Release Tool sheath is fully advanced.
- Do not re-sterilize or reuse the device. Reuse of the ReActiv8 Release Tool may result in elevated rates of the following risks: tissue damage, infection, and extended surgical procedure times.
- Do not use if the sterile package or components show signs of damage, if the sterile seal is ruptured, or if contamination is suspected for any reason.
- Do not use the tool if the expiration or "use-by" date has expired.
- This product contains up to 0.13% Cobalt (CAS No. 7440-48-4; EC No. 231-158-0). Prolonged exposure to metallic Cobalt is known to cause carcinogenic, mutagenic, and reproductive harm.
- In addition to risks commonly associated with surgery, the use of the ReActiv8 Release Tool includes, but is not limited to, the following risks: allergic reaction, tissue damage, infection, surgical complications, lead fragmentation, and oncogenesis.

Storage & Handling

- The Release Tool is a sterilized single-use device. Sterilization is completed via radiation.
- The Release Tool should be stored between -30°C to 60°C (-22°F to 140°F).
- Dispose of the device according to local guidelines/regulations.
- Only trained personnel may use this device.

Directions for Use

The procedure described below is recommended to remove implanted peripheral lead components. Fluoroscopic imaging capabilities may facilitate the procedure.

1. Disconnect leads from system components and dissect any fibrotic tissue to expose the proximal body of the lead. Perform any small fascial incisions as required to assist during dilation. *Note: A lead stylet is optional and may be used for added stiffness or lead handling.*
2. Advance the sheath down the lead, maintaining alignment, up to beginning of the proximal tine. A twisting motion with minimal axial force may assist sheath advancement past any proximal lead electrodes. The sheath will stop at the proximal tine.
3. Insert the dilator over the lead and fully inserted sheath and advance it up to and over as much of the proximal tine(s) as possible. **Note:** Take care not to use excessive axial force or push dilator and lead into other tissues, instead carefully using a twisting motion to advance the dilator along the lead and sheath.

4. Insert the trephine over the dilator, sheath, and lead. Use careful twisting motions to advance the trephine cutter while maintaining position of the dilator. A depth mark on the dilator will be seen when the trephine tip is aligned at the end of the dilator.
5. Use a gentle reciprocal rotation of the trephine with minimal axial load to incise any encapsulating tissue. Maintain axial alignment with the lead, periodically giving the lead a gentle tug to determine if it has been released. If released, the lead, dilator, and sheath should all pull out of the trephine cutter.
6. Once the leads have been removed, inspect them carefully to confirm that all portions of the lead and any tines have been removed.

Information and Support

Any serious incident that has occurred in relation to the device should be reported to the manufacturer below and the competent authority of the Member State in which the user and/or patient is established.

If there are questions regarding the use of the ReActiv8 Release Tool, please contact:

Mainstay Medical Limited
Clonmel House
Forster Way
Swords, Co. Dublin, K67F2K3,
Ireland

United States:
MML US, Inc.
2159 India Street,
Suite 200
San Diego, CA 92121

UK Responsible Person:
ICON (LR) Limited
100 Longwater Avenue
Green Park, Reading RG2 6GP
Berkshire, United Kingdom

Internet: www.mainstaymedical.com
E-mail: info@mainstaymedical.com
Tel.: (877) 702-8488

For label symbol definitions, refer to: www.mainstaymedical.com/symbols

Mainstay Medical and ReActiv8 are registered marks of Mainstay Medical Limited.

This product, and methods of use thereof, is covered by one or more of the patents identified at www.mainstaymedical.com/patents. This document and statement serve as notice under 35 U.S.C. § 287(a) with respect to patent marking.

Copyright © 2025 by Mainstay Medical Limited. All Rights Reserved. No portion of this IFU or manual may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording, or any information storage and retrieval systems, without the express written permission of Mainstay Medical Limited.

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Copyright © 2025 som tillhör Mainstay Medical Limited. Med ensamrätt. Ingen del av denna bruksanvisning får kopieras eller överföras i någon form eller på något sätt, elektroniskt eller mekaniskt, inklusive fotokopiering, inspelning eller via annat system för informationslagring och dataåtervinning, utan uttrycklig och skriftlig tillåtelse från Mainstay Medical Limited.
Copyright © 2025 - Mainstay Medical Limited. Tutti i diritti riservati. Nessuna parte di questo manuale può essere riprodotta o trasmessa in qualsiasi forma o con qualsiasi mezzo, elettronico o meccanico, inclusa fotocopiatura, registrazione o sistemi di registrazione e recupero delle informazioni, senza l'esplicita autorizzazione scritta di Mainstay Medical Limited.
MML VS, Inc.
Mainstay Medical en ReActiv8 zijn geregistreerde handelsmerken van Mainstay Medical Limited.