

Document downloaded from: <https://www.mainstaymedical.com>

Document Information:

Device Name: ReActiv8 Release Tool

Document/Label Type: Symbol Glossary









Country: US











Language: English



Revision: B

Filename: ReActiv8® Release Tool-Symbol Glossary-US-English- Rev B.pdf

ReActiv8 Release Tool Symbol Glossary

Symbol	Title and Reference Number	Title and designation number of standard	Definition
	Manufacturer 5.1.1	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the medical device manufacturer
	Medical Device 5.7.7	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the item is a medical device
	Catalogue Number: 5.1.6	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the manufacturer's catalogue number so that the medical device can be identified
	Batch code 5.1.5	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the manufacturer's batch code so that the batch or lot can be identified
	Unique device identifier 5.7.10	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates a carrier that contains unique device identifier information
	Date of Manufacture 5.1.3	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the date when the medical device was manufactured (yyyy-mm-dd)
	Use-by date 5.1.4	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the date after which the medical device is not to be used. Expiration Date (yyyy-mm-dd)
	Single sterile barrier system with protective packaging inside 5.2.13	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates a single sterile barrier system with protective packaging inside

	Do not re-use 5.4.2	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates a medical device that is intended for one single use only or for use on a single patient during a single procedure
	Sterilized using irradiation 5.2.4	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates a medical device that has been sterilized using irradiation
	Consult instructions for use or consult electronic instructions for use 5.4.3	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the need for the user to consult the instructions for use
	Do not use if package is damaged and consult instructions for use 5.2.8	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates that a medical device that should not be used if the package has been damaged or opened and that the user should consult the instructions for use for additional information
	Temperature limit 5.3.7	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the temperature limits to which the medical device can be safely exposed
	Humidity Limitation 5.3.8	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates the range of humidity to which the medical device can be safely exposed
	Keep dry 5.3.4	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates a medical device that needs to be protected from moisture
	Keep away from sunlight 5.3.2	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates a medical device that needs protection from light sources
	Contains hazardous substances 5.4.10	ISO 15223-1:2021(E): Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Indicates a medical device that contains substances that can be carcinogenic, mutagenic, reprotoxic (CMR), or substances with endocrine disrupting properties
	Protective Packaging is Recyclable U+2672	Unicode Standard ISO/IEC 10646:2020	Used as generic symbol for recycling or to indicate that material is recyclable
Rx Only	Prescription Use	N/A	Use by prescription only

	CE Mark	N/A	N/A
	UKCA Mark	N/A	N/A